



Better Science, Better Ethics: Animal-Free Cosmetic Safety Testing Approaches

RSC FTSG Personal Care Formulating
17 March 2022

Dr Louise Suleman
Project Manager, XCellR8



Coming up...

1. XCellR8 and our mission
2. Truly animal-free testing – what is it and why does it matter?
3. Adapting a test to animal-free conditions (TG442d)
4. Developing an animal-free test (XtraMild)





Our talented team

- UK-based GLP accredited lab
- Expertise in *in vitro* testing for safety & efficacy
- Regulatory and non-regulatory testing
- Fully customisable methods
- The only 100% animal-product-free laboratory globally



Our mission

To accelerate the world's transition to **100% animal-free testing** through our scientifically advanced and ethical approach



Are cosmetics now tested without cruelty to animals?

- Ban in Europe fully in place since 11th March 2013
- 80% countries still test cosmetics on animals!
- Other regulations (e.g. REACH) require animal testing of ingredients in some cases.
 - Recent rulings by ECHA have put the ban at risk in the EU.
- Many “non-animal” tests still ultimately require animal sacrifice
- Ecotoxicology testing still requires animals, e.g. acute toxicity to fish






Truly animal-free testing

- Most *in vitro* methods use animal components
 - Foetal bovine serum
 - Tissue extracts
 - Antibodies
- Scientific and ethical considerations
- Improved reproducibility (when using chemically defined systems)
- Driven by consumer and industry demand for sustainable, ethical products (*and* ethical testing)
- Truly animal-free testing needs to be animal-*product*-free
- Vegan products require vegan testing





The XCellR8 scale for animal-free testing



©XCellR8

Adaptation of an existing method:

KeratinoSens™

Skin Sensitisation Test (OECD TG 442d)

to Xeno-Free Conditions



Skin sensitisation – why testing is so important

- 1-3% (~7.5 million people) Europeans suffer from contact allergy to a cosmetic ingredient
- Skin sensitization to a cosmetic ingredient is a **permanent condition**
- *In vitro* tests provide an ethical alternative to human trials
- Preservatives and fragrances are the most common causes of skin sensitization in cosmetic products
- This is true for both natural and synthetic ingredients - **natural does not mean safe!**



Animal-product-free skin sensitisation testing


THE CHALLENGE

Current regulatory
guidance favours
“2 out of 3” approach

- DPRA (OECD TG 442c)
- KeratinoSens™ (OECD TG 442d)
uses animal components
- h-CLAT (OECD TG 442e)
uses animal components



Skin sensitisation adverse outcome pathway (AOP)



KEY EVENTS IN SKIN SENSITISATION AND RELATED TESTS

1. Contact
(Direct Peptide Reactivity Assay – **DPRA**)
2. Release of Pro-Inflammatory Cytokines by Keratinocytes (**KeratiNoSens™**)
3. Dendritic Cell Activation/Maturation
(human Cell Line Activation Test – **h-CLAT**)
4. Migration
5. T-cell Proliferation
(Local Lymph Node Assay - LLNA)

© XCellR8 Ltd

Regulatory guidance: “2 out of 3” approach

Adaptation of the KeratinoSens™ Skin Sensitisation Test (OECD TG 442d) to Xeno-Free Conditions

Published in ALTEX:



Belot, N., Sim, B., Longmore, CL., Roscoe, L. and Treasure, C. (2017)

[Adaptation of the KeratinoSens™ skin sensitisation test to animal-product-free cell culture >](#)



KeratinoSens™ - Method outline

- Human keratinocyte cell line (HaCaT) transfected with a luciferase reporter linked to Nrf2-mediated activation of Antioxidant Response Element (ARE)-linked genes
- 12 concentrations of test chemical incubated for 48 hours (in triplicate; 3 independent runs)
- Luciferase response measured by luminescence and cytotoxicity measured by MTT





Xeno-Free adaptation of KeratinoSens™

- Animal-derived components were replaced with human-derived & recombinant equivalents:
 - FBS replaced with pooled **human serum** (60-70 donors) obtained from FDA-approved source / Sigma Aldrich – cells adapted to new culture conditions
 - Porcine trypsin replaced with **recombinant Trypzean™**
- In-house validation using the panel of proficiency chemicals and performance standards for OECD TG 442d



Results: Non-Sensitisers (*as per LLNA*)

Chemical Name	Validated Reference Method (VRM)			XCellR8 Animal-Product-Free Adaptation		
	I _{Max}	EC1.5 (µM)	Prediction	I _{Max}	EC1.5 (µM)	Prediction
Isopropanol	1.2	n.i.	Non-Sensitiser	1.2	n.i.	Non-Sensitiser
Salicylic Acid	1.1	n.i.	Non-Sensitiser	1.4	n.i.	Non-Sensitiser
Lactic Acid	1.3	n.i.	Non-Sensitiser	1.3	n.i.	Non-Sensitiser
Glycerol	1.2	n.i.	Non-Sensitiser	1.4	n.i.	Non-Sensitiser
<u>4-methoxy-acetophenone</u>	1.7	449.3	<i>Sensitiser</i>	2.1	620	<i>Sensitiser</i>
Chlorobenzene	1.2	n.i.	Non-Sensitiser	1.2	n.i.	Non-Sensitiser
Methyl Salicylate	1.2	n.i.	Non-Sensitiser	1.2	n.i.	Non-Sensitiser
Sulfanilamide	1.4	n.i.	Non-Sensitiser	1.1	n.i.	Non-Sensitiser

© XCellR8 Ltd

n.i. = not induced



Results: Sensitisers (as per LLNA)

Chemical Name	Validated Reference Method (VRM)			XCellIR8 Animal-Product-Free Adaptation		
	I _{Max}	EC1.5 (µM)	Prediction	I _{Max}	EC1.5 (µM)	Prediction
Cinnamyl alcohol	1.7	123.6	Sensitiser	4.2	20	Sensitiser
Ethylene Glycol Dimethacrylate	188	57.4	Sensitiser	4.8	29	Sensitiser
<u>Phenyl Benzoate</u>	1.3	n.i.	<i>Non-Sensitiser</i>	1.1	n.i.	<i>Non-Sensitiser</i>
<u>Eugenol</u>	1.3	n.i.	<i>Non-Sensitiser</i>	2.2	286	<i>Non-Sensitiser (borderline)</i>
2-Mercaptobenzothiazole	8.8	48.1	Sensitiser	6.9	57	Sensitiser
Citral	96.4	23.2	Sensitiser	3.8	18	Sensitiser
Isoeugenol	6.4	16.1	Sensitiser	3.4	20	Sensitiser
Methyldibromo Glutaronitrile	4	7.8	Sensitiser	2.7	8	Sensitiser
4-Methylaminophenol Sulfate	5.9	9.4	Sensitiser	36.1	4	Sensitiser
Para-phenylene Diamine	26.8	5	Sensitiser	28.2	6	Sensitiser
2,4-Dinitrochlorobenzene	14.8	2.5	Sensitiser	8.5	1	Sensitiser
4-Nitrobenzyl Bromide	6.9	1.3	Sensitiser	10.5	<0.98	Sensitiser
Oxazolone	2.4	175.5	Sensitiser	5.4	129	Sensitiser

n.i. = not induced



Animal-product-free (APF) adaptation of KeratinoSens™ Conclusions

- All 20 reference chemicals correctly classified in line with Validated Reference Method (VRM)
- Data accepted by the OECD Expert Working Group on Skin Sensitisation and WNT National Co-Ordinators' Committee
- Adapted method published as an Annex to the VRM in the new version of OECD TG 442d 2018
- Therefore full acceptance as a regulatory method






- Participation in thought-starter paper and OECD workshop on the ethical use of human reagents:
 - [Addressing potential ethical issues regarding the supply of human-derived products or reagents in *in vitro* OECD Test Guidelines.](#)

Published in ALTEX 2019

- Xeno-free adaptation of h-CLAT including human serum and animal-free antibodies:
 - [Edwards *et al*, ALTEX, 2018](#)
- Adaptation of KeratinoSens™ and h-CLAT to fully defined conditions



Creating an animal-free test:

XtraMild – Skin Mildness Test for Safety &
Claim Support



Why we need a new method to predict mildness (I)

- Study of 12,377 individuals in Europe*
- Incidence of skin reactions lasting more than 3 days:
 - 19.3% within the last month
 - 31.8% within the last year
 - 51.7% within a lifetime
- Avoidance of daily life consumer products due to skin reactions:
 - 37.0% for skincare
 - 17.7% for “household or functional” products

* Naldi *et al* (2014). Prevalence of self-reported skin complaints and avoidance of common daily life consumer products in selected European regions. *JAMA Dermatol* **150**(2): 154-162





Why we need a new method to predict mildness (II)

- Increasing demand from consumers for ever milder products, that they feel confident using even when their skin is feeling extra sensitive
- Increasing demand from marketing teams for differentiating claims
- 2 year research project 2017-2019, funded by Innovate UK
- Research aims:
 - Optimising *in vitro* and human *in vivo* test methods for maximum sensitivity
 - Assess predictive capacity
 - Real world applications


We work with
Innovate UK





Existing methods: *In vitro* irritation testing


- 3D human skin models, grown at the air-liquid interface
- Suitable for testing ingredients and finished products
- Applied directly to the tissue surface – good model of “real life” exposure
- Standard regulatory method (OECD TG 439) measures a single exposure time to classify irritants vs non-irritants for hazard identification and labelling purposes
- Validated against historical animal data (Draize test)
- A more sensitive approach is required for today’s mild cosmetic ingredients and formulations beyond a yes/no answer – how mild is the test item?



Cross section through reconstructed human epidermis



The ET50 method



- Measures cell damage over a time course
- Classifies as Severe, Moderate, Mild or Minimal / Non-Irritant
- ET50 = time taken to reduce the viability of the skin model to 50% compared with untreated controls
- ET50 values allow rank order of irritation to be determined in comparison with other formulations / competitor and market leading products
- Standard methodology limited to 18 hours




How we optimised the test methods *in vitro*

- Development of an extended timepoint *in vitro* 3D model to look at the irritancy potential of ultra-mild test items over 48 hours
- Determination of ET₅₀ values for known surfactant controls with a range of irritation potentials
- Development of a prediction model linking the *in vitro* skin irritation ET₅₀ method with an *in vivo* human skin patch test model for ultra-mild surfactants
- Creation of a database of industry leading ingredients and formulations to be used as benchmarks in future tests for client companies



How we optimised the test methods *in vitro*

Test Items	Surfactants: SLS, SLES, CAPB, a novel “mild” surfactant Applied to the skin model surface and incubated for 1, 5, 18, 24 and 48 hours
Controls	Negative control: not treated Positive control: Triton X-100 (non-ionic surfactant): 1% solution
Measurement	Metabolic activity (conversion of MTT) as an indicator of cell damage
Output	ET50 value (time taken to reduce the viability of the cells to 50% compared with the untreated negative control)



Tissue Culture insert


Determining the correlation between *in vitro* and *in vivo* results – some examples

1. SURFACTANTS
2. FACE MASKS




In vitro irritation potential of 4 surfactants

ET₅₀ determination of A
0.3% SLES



ET₅₀ determination of C
0.3% SLS




Test items (0.3%, pH 4.7) at 1, 5, 18, 24, 48hrs


Rank order of irritancy using linear extrapolation and logic equation

	C	>	A	>	B	>	D
ET50	9.37		10.25		29.4		38.08
CS	14		9		4		0

ET₅₀ determination of B
0.3% CAPB



ET₅₀ determination of D
0.3% Novel surfactant




Irritancy classification:

- C = SLS:** Moderate to Mild
- A = SLES:** Moderate to Mild
- B = CAPB:** Non-Irritant
- D = Novel surfactant:** Non-Irritant



Face mask comparison *in vitro* and *in vivo*

ET₅₀ determination of 3 face mask formulations



Rank order of irritancy using linear extrapolation and logic equation

	B	>	A	>	C
ET50	12.86		14.42		>48
CS	11		5		2

IRRITANCY CLASSIFICATION


B = face mask 2: Very mild

A = face mask 1: Very mild

C = face mask 3: Non-irritating

Face mask C is the mildest product using this method

Building a reference database: relative mildness of soaps vs facial cleansers





A variety of applications

- Ingredients:
 - Assessment of novel biosurfactants and other ingredients to assess mildness compared with other manufacturers and traditional materials
- Formulations:
 - *In vitro* benchmarking of new products against other brands or in-house formulations in development
 - Growing database for benchmark values currently includes:
 - ✓ Facial soaps
 - ✓ Facial cleansers
 - ✓ Face masks
 - ✓ Moisturisers
 - ✓ Body soaps
 - ✓ Shampoos
 - ✓ Shower gels
 - ✓ Sunscreens
 - ✓ Deodorants
 - ✓ Baby care products (oils, lotions, shampoos, bubble baths)





Take home...

- Many 'non-animal' tests still use animal components
- Using animal components raises scientific & ethical questions
- Adaptation of regulatory safety tests is already happening but it needs investment
- Technology and expertise to develop new animal product free methods is available
- *In vitro* testing provides robust safety information, and strong database for bench marking ingredients & formulations – IT'S BETTER SCIENCE
- Ethical advantages: limits human exposure, whether used as stand-alone test or pre-screen to clinical studies
- Marketing / consumer appeal

Thank you

Dr Louise Suleman
Project Manager
louise.suleman@x-celr8.com

